

## Product Specification

Ref: **4891-531**

Alias: 489199531

Brand:

**ZILLION<sup>®</sup>**  
**BLACK**

DRQ 14/21 REV.02 2024/09/20

<b>Description:</b>	Sponge forceps
<b>Final dimensions:</b>	26.5 cm
<b>Composition:</b>	High performance polymer.
<b>Latex:</b>	Free
<b>PVC:</b>	Free
<b>Re use:</b>	No
<b>CMR:</b>	Free
<b>Endocrine Disruptors:</b>	Free
<b>Colour:</b>	Grey
<b>Sterilization method:</b>	Ethylene Oxide
<b>Shelf life:</b>	60 months
<b>CFU/g</b>	<10 <sup>2</sup>

STERILE EO






isento  
Latex  
free

isento  
PVC  
free



<b>Nr and Revision Date:</b>	44 -- 2024/02/28
<b>CE Class:</b>	Ila
<b>CE mark:</b>	CE1639
<b>Compliance with:</b>	MDR (EU) 2017/745

<b>Applicable standards:</b>	EN ISO 11135:2014 ; EN ISO 20417:2021 ; EN ISO 14971:2019 ; EN ISO 15223-1:2021 ; ISTA Standard 2A ; ASTM D4169-22 ; EN ISO 10993-18:2020 ; EN ISO 10993-1:2020 ; EN ISO 11607-1:2020 ; EN ISO 11607-2:2020 ; EN ISO 14971:2019 A11:2021 ; ISO 2759:2014 ; ISO 3037:2013 ; EN ISO 11135:14 A1:2019 ; CEN ISO TS 16775 2021 ; ISO TR 24971:2020 ; UDI FAQs ; EN ISO 13485:2016 ; EN ISO 13485:2016 AC:2018 ; EN 556-1:2001 AC:2006 ; ISO 14644-1:2015 ; ISO 14644-5:2004 ; ISO 14698-1:2003 ; EN ISO 13485:2016 A11:2021 ; EN ISO 10993-7:2008 ; ISO 10993-18:2020 AMD1:2022 ; EN ISO 11737-1:2018 ; ISO 10993-7:2008 Amd1 2019 ; EN ISO 11737-1:18 A1:21 ; ISO 11737-3:2023 ; EN ISO 11737-2:2020 ; EN 17141:2020 ; ISO 14644-2:2015 ; ISO 14644-3:2019 ; ISO 14698-2:2003 ; ASTM D7386-16 ; ASTM F1980-21 ; ISO 8601-1:2019 ; ASTM F1886/F1886M-16 ; ASTM F1929-15 ; ASTM F88/F88M-21 ; IEC 62366-1:2015 A1:2020 ; IEC TR 62366-2:2016 ; ISO TR 20416:2020 ; ISO 2859 1:1999 ;
<b>GMDN:</b>	62490 - Dressing/utility forceps, scissors-like, single use

PACKING					
Level	Units	Type	Dimension	Stamp	Labels
	1	1 Rigid peel-open pack	200x320mm		
	2	9 Folding box	374x204x190mm		labeled
	3	81 Carton box	586x388x641mm		labeled

<b>Basic UDI</b>	56056221207002JV
	EAN 1th level/nivel:05605622057842
	EAN 2th level/nivel:15605622057849
	EAN Carton Box/Caixa:25605622057846

<b>Important Notes:</b>	The device is intended to be used on a single patient during one single procedure to prevent risks from contamination, including possible infections, and/or functional damage.
-------------------------	---

<b>Storage instructions</b>	Store in a dry place at room temperature and protect from sunlight. The product has a good stability and therefore, normal conditions of medical devices storage shall apply
-----------------------------	--

<b>Waste instructions</b>	Dispose according to local laws rules, regulations and/or infection prevention standards for healthcare waste, to minimize potential risk of contamination, infection or pollution.
---------------------------	---

<b>Intended purpose</b>	For holding swabs for absorption during surgery or woundcare procedures or for skin disinfection procedures.
-------------------------	--

*Any modification of the product as presented or any use out of its intended purpose are not of the manufacturer responsibility*

**ISSUED BY: Quality Management Department, 2024/11/21**

---

CONFIDENTIAL: The information contained herein is confidential and belongs to Bastos Viegas, s.a. It is provided for specific purposes and in connection with established agreements. It is not permitted disclosure, reproduction, copying or distribution, in whole or in part, for other purposes or to any third parties.

---