

Product Specification

Ref: 4891-501

Alias: 489199501

Brand:

ZILLION[®]
BLACK

DRQ 14/21 REV.02 2024/09/20

Description: Curette


Composition:	High performance polymer.
Latex:	Free
PVC:	Free
Re use:	No
CMR:	Free
Endocrine Disruptors:	Free
Colour:	Grey
Sterilization method:	Ethylene Oxide
Shelf life:	60 months
CFU/g	N.A.
Sterility Assur. Level_SAL	≤10 ⁻⁶
Length	139,5 mm
Break test	>8 N

STERILE EO



isento
Latex
free

isento
PVC
free



Nr and Revision Date:	48 -- 2024/08/30
CE Class:	Ila
CE mark:	CE1639
Compliance with:	MDR (EU) 2017/745

Applicable standards:	EN ISO 10993-18:2020 ; EN ISO 10993-1:2020 ; ISTA Standard 2A ; ASTM D4169-22 ; EN ISO 15223-1:2021 ; EN ISO 20417:2021 ; EN ISO 14971:2019 ; ISO 3037:2013 ; EN ISO 11607-1:2020 ; EN ISO 11607-2:2020 ; CEN ISO TS 16775 2021 ; EN ISO 11135:14 A1:2019-Not applicable ; EN ISO 14971:2019 A11:2021 ; EN ISO 13485:2016 A11:2021 ; ISO 2759:2014 ; ISO TR 24971:2020 ; UDI FAQs ; EN ISO 13485:2016 ; EN ISO 13485:2016 AC:2018 ; EN 556-1:2001 AC:2006 ; ISO 11135:2014 ; ISO TS 16775:2014 ; ISO 14644-1:2015 ; ISO 14644-5:2004 ; ISO 14698-1:2003 ; EN ISO 10993-7:2008 ; ISO 10993-7:2008 Amd1 2019 ; ISO 11737-3:2023 ; ISO 10993-18:2020 AMD1:2022 ; EN ISO 11737-2:2020 ; EN ISO 11135:2014 ; EN 17141:2020 ; ISO 14698-2:2003 ; ISO 14644-2:2015 ; ISO 14644-3:2019 ; ASTM D7386-16 ; ASTM F1886/F1886M-16 ; ASTM F1929-15 ; ASTM F88/F88M-21 ; ASTM F1980-21 ; ISO TR 20416:2020 ; ISO 2859 1:1999 ; ISO 8601-1:2019 ; IEC 62366-1:2015 A1:2020 ; IEC TR 62366-2:2016 ;
GMDN:	62558 - General-purpose curette

Important Notes: Rigid peel-open pack: made of Tyvek and rigid film

PACKING

Level	Units	Type	Dimension	Stamp	Labels
4	1	1 Rigid peel-open pack	53x173mm		
20	2	30 Folding box	176x056x212mm		labeled
20	3	120 Carton box	245x188x223mm		labeled

Basic UDI	56056221202002HS
	EAN 1th level/nivel:05605622241258
	EAN 2th level/nivel:05605622441252
	EAN Carton Box/Caixa:05605622041254

Important Notes: The device is intended to be used on a single patient during one single procedure to prevent risks from contamination, including possible infections, and/or functional damage.

Storage instructions Store in a dry place at room temperature and protect from sunlight.The product has a good stability and therefore, normal conditions of medical devices storage shall apply.

Waste instructions Dispose according to local laws rules, regulations and/or infection prevention standards for healthcare waste, to minimize potential risk of contamination, infection or pollution.

Intended purpose For scraping, debriding and/or cleaning skin and soft tissues and wounds.

Any modification of the product as presented or any use out of its intended purpose are not of the manufacturer responsibility

Environment - Packaging waste	Weight (gr) of packing material used in 100 pcs	981
-------------------------------	---	-----

ISSUED BY: Quality Management Department, 2024/11/21

CONFIDENTIAL: The information contained herein is confidential and belongs to Bastos Viegas, s.a. It is provided for specific purposes and in connection with established agreements. It is not permitted disclosure, reproduction, copying or distribution, in whole or in part, for other purposes or to any third parties.